



Wellmed Pro | Wellview Universal Group

TECHNICAL DATA SHEET

Automatic Erythrocyte Sedimentation Rate Analyzer

YSESR-2010 autotiered tidal titration analyzer

Technical Parameter

Analysis Time

30 min

Analysis Capacity

10 tests/30min

Reading Channel

10

Load Capacity

Max of 10 samples/batch

Load Form

Measure at any time

Analysis Result

Westergren results (mm/h)

Measurement Methods

Infrared detection, automatic temperature compensation and time correction

Reading Resolution

0.2 mm

Result Resolution

1 mm

Range of Liquid Level

50 mm - 64 mm

Display

LCD display

Operation

Touch screen

Interface

RS232 serial interface

Software

Complies with ISO15189 standard

Coefficient of variation

less than 15%

Storage

No less than 500 specimens

Protection

Power-off protection

Print

External. Reports with the curves can be printed while connecting to the computer (nonstandard configuration).

Size (height × length × width)

300 mm × 400 mm × 300 mm

Weight

5.65 kg

Power supply

AC 220 ± 22 V, 50 Hz ± 1 Hz

Working conditions

Temperature: 15 C - 32 C; relative humidity: 45% - 80%

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.