



Wellmed Pro | Wellview Universal Group

TECHNICAL DATA SHEET

Advanced Insulin Injection Training Pad Model YSBIX-HS17

Cloth elastic can be worn in the specific site of the patient for self-injection of insulin drug, such as: abdomen, thigh, upper arm, etc. The length and tightness can be adjusted. And the thickness of the pad is suitable for different specifications of injection syringe. The hard plate in the back would prevent wearing situation of pad.

Features:

- 1.Imported plastic material, the skin soft and flexible, with a cloth belt, can be worn on the body of a patient specific parts, such as the abdomen, thigh, arm etc..
- 2.Insulin drug self injection. Module thickness allows the use of different specifications of the puncture syringe.
3. For security reasons, on the back of this module ovoid with hard plastic board, prevent injection practice module stabbing injury to practice.
4. The elastic band length and tightness can be adjusted

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.