



Wellmed Pro | Wellview Universal Group

TECHNICAL DATA SHEET

Autoclavable Multi-Channel Laboratory Pipette YSTE-YY2P

- Dispensing head rotates for effortless pipetting convenience
- Individual piston and tip cone assembly allow for easy repair and maintenance
- Compound material-made tip cone secures high sealing performance
- Compatible with most universal tip brands

Volume Range

Increment

Test Volume(μ l)

Accuracy error

Failed accuracy

%

μ l

%

μ l

0.5-10 μ l

0.1 μ l

10

1.50%

0.15

1.50%

0.15

5

2.50%

0.125

2.50%

0.125

1

4.00%

0.04

4.00%

0.04

5-50 μ l

0.5μl
50
1.00%
0.5
0.50%
0.25
25
1.50%
0.375
1.00%
0.25
5
3.00%
0.15
2.00%
0.1
50-300μl
5μl
300
0.70%
2.1
0.25%
0.75
150
1.00%
1.5
0.50%
0.75
50
1.50%
0.75
0.80%
0.4
Volume Range
Increment
Test Volume(μl)
Accuracy error
Failed accuracy
%
μl

%
μl
0.5-10μl
0.1μl
10
1.50%
0.15
1.50%
0.15
5
2.50%
0.125
2.50%
0.125
1
4.00%
0.04
4.00%
0.04
5-50μl
0.5μl
50
1.00%
0.5
0.50%
0.25
25
1.50%
0.375
1.00%
0.25
5
3.00%
0.15
2.00%
0.1
50-300μl
5μl
300
0.70%

2.1
0.25%
0.75
150
1.00%
1.5
0.50%
0.75
50
1.50%
0.75
0.80%
0.4

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.