



Wellmed Pro | Wellview Universal Group

TECHNICAL DATA SHEET

Automated Tissue Processing Machine YSPD-TS150

- © English menu, man-machine interaction, simple operation.
- © Prompt judge and automatic processing in case of abnormal condition, which can prevent the tissue from being damaged.
- © Double-control heat without water
- © Based on different tissues, several running programs available for customers
- © Self-provided power supply can ensure to continue 6 hours' working even if the electricity has been cut off.
- © Environment friendly: The air as purified through internal circulation can be exhausted off the window.
- © With transmission of gear and rack, mechanical failures can be reduced.
- © Several standby programs available

Technical Parameter:

- © 12 cylinders (3 wax cylinders in the end)
- © Capacity of single-cylinder: 1500ml
- © Basket can hold 90 specimens (Standard dehydrating box)
- © Single-cylinder processing time: Randomly set within the range from 0 to 99 hour(s).
- © The temperature of wax cylinder can be randomly set within the range of 99°C.
- © Time for cylinder to be lifted: 38 seconds
- © Tissue-protection cylinder sequence: 7th cylinder
- © Times of mixing per minute: 3 times

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.