

TECHNICAL DATA SHEET

CHISON ECO1 Ultrasound Machine



Reference photograph. Confirm accessories and configuration in the quotation.

CHISON ECO1 Ultrasound Machine For Sale With Best Price

User-friendly and Modern Design

Support PC Printer

Advanced Technologies

SRA is the technique that uses a variety of denoising algorithm to suppress speckle, smooth the tissue images, and

make the faint edge more appearing.

Applications

Abdominal , OB/GYN , Urology, Pediatric; Small Organ (breast, testes, thyroid);

Cardiac, Vascular Peripheral, Musculo-Skeletal Conventional & Surface

Display Modes

B, B/B, 4B, M, B/M

Zoom

PIP zoom & pan zoom

Advanced Technologies

THI, i-Image, SRA, Compound

Option

2nd Probe connector , DICOM 3.0, Water-proof keyboard cover, Carrying Case (BG-100) , Support PC Printer, TR-9000 , Built-in battery

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.