

TECHNICAL DATA SHEET

Electrosurgery Generator YSESU-D5



Reference photograph. Confirm accessories and configuration in the quotation.

Electrosurgery Generator

Tissue Monitoring Technology

This technology monitors the tissue density, especially for polypectomy & papillotomy, cutting is smartly divided incut & pause, so that coagulation occurs in each pause.

The system is intelligent to change the pause interval as tissue density changes, hence produces fine cutting with coagulation effect.

Patient Return Electrode Monitoring (PREM)

REM (Return Electrode Monitoring) function monitors the range of contact area between a patient and the neutral plate. If the size is inappropriate, it automatically stop the high-frequency current with audible alarm to minimize the danger of burning incidents.

Compatible with other surgical equipment

Endoscopy, TURP, HV Argon Plasma Coagulator, Optimum Smoke Evacuation System.

Waterproof Panel Operation

Microcomputer Auto Control, Three-zone (LED) digital display.

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.