

TECHNICAL DATA SHEET

Laboratory Roller Mixer Price YSTE-MR2A

YSTE-MR2A Laboratory Roller Mixer

Features:

- ⊗ It can carry out two functions(Roll and Swing), adapted to mixing of blood samples, adhesive materials and Liquid - solid suspensions.
- ⊗ Driven by High performance DC declaration Motor, Running stable, No noise, Long lifetime.
- ⊗ 7 rolls, compact structure makes it less space occupation. adapted for obvious environment.
- ⊗ LED with backlight can display the current speed, time and operation mode, convenient for reading and setting.
- ⊗ Precise digital time and speed setting can easily carry out the continuous operation mode.
- ⊗ Wide-range voltage design. Avoid the interference to the equipment by unstable voltage .
- ⊗ Reasonable space between the rolls makes it fit for different size tubes.
- ⊗ According to actual requirements, we can supply corotation-reverse cross model, The roller can be lengthened within a certain range to fit for more tubes.

Main Technical Data:

Operation method:

Swing angle:

Rating speed:

Max. Loading weight:

Running mode:

Timing range:

Input voltage(frequency):

Permissible relative moisture:

Output power:

Roll length:

Protective class:

Permissible ambient temperature:

Dimensions:

Net. Weight:

Roll+swing

22mm±1mm

20-100rpm

1.5kg

Constant, inching
0-999min
AC100-240V(50/60Hz)
80%
10W
236mm
IP21
5-50°C
410×270×100mm
4.0kg

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.