



Wellmed Pro | Wellview Universal Group

TECHNICAL DATA SHEET

18L Benchtop Class B Dental Autoclave

18L Benchtop Class B Dental Autoclave

Product parameters:

Model

Voltage/Power

Chamber Size

(mm)

Outside Size

(mm)

Packing Size

(mm)

Net Weight

(kg)

YSMJ-TZO-E18

220V/1500W

φ245×350

471×586×395

660×525×495

47.4

YSMJ-TZO-E23

220V/1700W

φ245×450

471×696×395

770×530×495

53.4

Product Function:

1. The main system specially revised according to EN13060 standard with intelligent.

Adjusting sterilizing time and vacuuming time, more accurate control for stable and crashing sterilization effect.

2. New LCD display window, 5 buttons control panel, for convenient and easy operation.

3. Intelligent help system could demonstrate the condition of sterilization automatically.

4. Both printer and USB output, you could choose the way which you like to backup the sterilization records.

The specification on the output is more accurate and detailed.

5. Memory system inside, you could check the sterilization.

6. Adding a new set temperature detecting device. There are two sets separate temperature detecting device inside.

One is for controlling the temperature, the other one is for checking the controlling temperature, it could prevent

if the temperature sensor for controlling temperature disabled or incorrect and bring the incompletely sterilization.

The precision of temperature control developing decuple, it could be precision to 0.1 degree.

7. Adjusting the check mode of vacuuming, adopt higher parse defree pressure sensor, the precision could be 0.001bar,

it will collect the temperature or pressure value per 30ns a time during the program. More accurate value collects for

exact analysis and self-checking automatically.

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.