



Wellmed Pro | Wellview Universal Group

TECHNICAL DATA SHEET

Manual door vacuum steam sterilizer XG1.C

XG1.C manual door vacuum steam stabilizer

Product introduction

It is applicable for the steriliation of unpacked solid instrument in department of stomatology, ophthalmology, operation room and CSSD.

Product features

.Touch screen display of temperature, pressure, time, running status and failure warning.

.Equipped with professional and systematized sterilizing programs of fabric, instrument, rubber and liquid.

.The positive pressure air-discharging system will effectively discharge the cold air from chamber to improve steam saturability.

.Build-in efficient steam generator, microcomputer temperature control, time and energy saving.

.Equipped with post-drying function which makes the residual humidity of sterilized articles less than 3%.

Equipped with printer to print date, time and process parameters.

.Over-heat auto protection device, dual door safety lock device, overpressure relief valve, short-circuit protection and overpressure protection.

.Perfect four-level authority management system to avoid misoperation.

632x1000mm

1400*890x1780mm 750kg AC 380V,50Hz

17Kw

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.