

TECHNICAL DATA SHEET

MINDRAY BC-5380 5 Part Diff CBC Auto Blood Hematology Analyzer

mindray



Reference photograph. Confirm accessories and configuration in the quotation.

Model	BC-5380
Parameters	27
Histograms / scattergrams	3 / 1
Sampling mode	Autoloader; closed tube
Autoloader throughput	Up to 60 samples/hour
Closed-tube throughput	Up to 50 samples/hour
Whole-blood sample volume	20 μ L
Prediluted sample volume	20 μ L
Test modes	CBC; CBC+DIFF
Data storage capacity	Up to 40,000 results including numeric and graphical information
Communication	LAN port supports HL7 protocol
Operating temperature	15~30 °C
Operating humidity	30~85%
Operating air pressure	70~106 kPa
Power requirement	AC 100-240 V; $\leq$ 300 VA; 50/60 Hz
Weight	$\leq$ 65 kg

Exact BC-5380 own technical sheet visually verified. Tiny performance table and analyte abbreviations omitted rather than relying on OCR.

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.