

TECHNICAL DATA SHEET

Slide Microtome YSPD-Q205



Reference photograph. Confirm accessories and configuration in the quotation.

Slide Microtome YSPD-Q205

- © It have clipper -built case , the user can put blade , wax block ,or other consumable one the top of case . it is easy .
- © it donot need to lubricate and maintain the microtome frequently .
- © The handwheel in the right can be locked in any position , and the handle can be locked in highest position , it is double safe for changing blade or tissue block .
- © The specimen Orienting system: axes X and Y 8°, axis 180°.
- © It adopt the special blade holder , it is more safe and convenient to use the disposable blade directly . it also can to longer the user life of disposable blade

Technical Data

- © Slice thickness range: 0- 60um
- Setting values : from 0-2 um in 0.5 um -increments
- Form 2-10um in 1 um-increments
- Form 10-20 um in 2 um-increments
- Form 20-60um in 5 um-increments
- © Horizontal especmen stroke: 28 mm
- © Vertical especmen stroke: 60mm
- © Slip precision: $\pm 5\%$
- © Maximum slice Section: 50 × 45mm

Panama | Central America | Caribbean | CARICOM

Information from the authorized catalog. Confirm configuration, accessories and documentation in the quotation.

Assisted translation: review the manufacturer specifications before purchasing.